

# Lab Notebook Entry 47

---

## Date Performed

2-Jun

## Procedure

miniprep pThermolyse backbone

## Details

Purpose: Get the pThermolyse backbone out of the cells in DNA form. Materials: NOTE: Our P1000 pipet doesn't work right now, so don't try using it, Miniprep kit, pThermolyse backbone starter cultures, eppendorf tubes

## Protocol

Do one miniprep for each of the starter cultures. Use the protocol for minipreps in the 2017 protocol file, EXCEPT: double the amounts of buffers P1, P2, and N3 used. Also preheat the buffer EB by sticking it in the incubator for a few minutes before you use it.