

Lab Notebook Entry 28

Date Performed

17-May

Procedure

Run new pTlpA PCR product on a gel, then do gel extraction

Details

Purpose: Try and get a higher concentration of pTlpA DNA than we did the first time.

Materials: Agarose, 6X purple loading dye, ladder, 1X TAE buffer, GelRed, pTlpA36 PCR product, empty PCR tube for each sample, Razor blade

Protocol

Use the protocol that is in the Gel Extraction kit. RUN THE WHOLE PCR PRODUCT ON THE GEL