

# Lab Notebook Entry 19

---

## Date Performed

5/9/17

## Procedure

Run a gel of 4th attempt of holin/endolysin PCR, 3rd attempt of pTlpA36 PCR, and 2nd attempt of antiholin

## Details

Purpose: Check that our PCR products are the correct size. Materials: Agarose, 6X purple loading dye, ladder, 1X TAE buffer, GelRed, pTlpA36-wasabi PCR product, 2nd attempt antiholin PCR product, 2nd attempt endolysin/holin PCR product, empty PCR tube for each sample

## Protocol

Use the protocol in the 2017 protocols file for the gel.